

# BorsoPleat-P

## Absolute Rated Pleated Polypropylene Cartridge Filters



**A range of absolute rated cartridge filters from Van Borselen Filters, featuring the latest developments in meltblown polypropylene filter media technology, BorsoPleat-P cartridges are based on a robust all polypropylene construction, offering removal ratings from 0.5 to 105 micron absolute.**

BorsoPleat-P cartridges are suitable for absolute removal of unwanted particulates and for prefiltration to membrane filter.

The graded multi-layer polypropylene media provide prefiltration of the process fluid prior to the absolute rated final layer. The unique design of the BorsoPleat-P cartridges helps to achieve lower running costs and a smaller process footprint.

The BorsoPleat-P are also highly resistant to integrity failure caused by steam sterilisation and have excellent chemical compatibility characteristics.

They are suitable for applications ranging from bioburden reduction and the clarification of a wide range of process liquids and end products.

## Applications

BorsoPleat-P cartridges provide absolute filtration where reproducibility and consistency of performance are critical. Suitable for the filtration of aqueous and organic liquids, BorsoPleat-P cartridges can be used as prefilters or final filters in the following applications:



- **Pharmaceuticals and Bioprocessing**

The batch preparation of intermediates used in the manufacture of pharmaceutical and bioprocessed products. Used as prefilters to extend the life of membrane filters.

- **Food and Beverages**

The clarification of beers, wines and spirits to a clear and bright finish without affecting taste or colour.

- **Inks and coating**

For removal of agglomerates and gels with minimal loss of product.

- **Fine chemicals**


The filtration of high grade chemicals including solvents, reagents, photographic emulsions and plating solution.

- **Cosmetics**

The clarification of process water and intermediates for the finished product.

- **Process Water Systems**

The filtration of process water installations for removal of general contamination and resin fines.



### Features and Benefits

#### • BorsoPleat-P cartridges

Extensive research and selection of the latest and most advanced polypropylene meltblown microfibre filter media, results in improved performance, leading to extended filter life at a given efficiency.

#### • Graded multi-layer media

The multi-layer media structure provides prefiltration of the process fluid prior to the absolute rated final layer. This combination provides economy of use and a smaller process footprint.

#### • High filtration area

Large surface area for low clean pressure drop.

#### • Guaranteed removal ratings

BorsoPleat-P cartridges are validated using the recognised industry standard modified OSU-F2 single pass test to Beta 5000 (99.98% efficiency).

#### • Excellent chemical compatibility

Resistant to many process chemicals, BorsoPes Biological cartridges are suitable for use in a wide range of process applications.

#### • Suitable for steam and hot water sanitisation

BorsoPleat-P cartridges are resistant to repeat steam sterilisation up to 135°C (275°F) and hot water cycles at up to 90°C (194°F).

#### • Environmentally friendly

BorsoPleat-P filters are environmentally friendly, all spent cartridges can be readily incinerated to trace ash.

#### • Full traceability

All BorsoPleat-P cartridges are identified with a batch serial number. Each BorsoPleat-P cartridge is supplied with a Certificate of Quality and an operating instruction leaflet.

#### • Controlled manufacturing environment

BorsoPleat-P cartridges are manufactured in an ISO Cleanroom environment by fully gowned staff, minimising the risk of contamination.

### Cartridge Construction

The high quality robust all polypropylene construction of BorsoPleat-P cartridges, allows for excellent chemical compatibility with a wide range of fluids.

The meltblown polypropylene media provides a bonded matrix thus eliminating fibre migration.

The inherent structural stability of the BorsoPleat-P, prevents 'channelling' and avoids the risk of particle unloading even under impulse conditions.


The multi-layer combination of filter media, irrigation mesh and drainage material carefully pleated and thermally bonded maximises the media area and ensures an efficient flow throughout the cartridge.

The BorsoPleat-P fusion bonded construction ensures cartridge integrity. No surfactants or bonding agents are used, minimising extractables.

Table 1 Particle Retention Rating

Code	Pore Rating (microns)	Absolute Rating 99.98% Beta 5000 (microns)	Nominal Rating 99% Beta 100 (microns)	Nominal Rating 90.00% Beta 10 (microns)
BLP	0.5	0.5	0.1*	<0.1
BLP	0.8	0.8	0.3	0.2*
BLP	1	1	0.4	0.2*
BLP	2	2	1.1	0.4
BLP	3	3	0.9	0.45
BLP	5	5	1.5	0.5
BLP	7	7	2	0.6
BLP	10	10	4	1
BLP	15	15	7.5	1.5
BLP	20	20	8	2
BLP	30	30	12.5	6
BLP	40	40	15	8
BLP	60	60	45	25
BLP	90	90	65	30
BLP	105	105	95	75

\* Extrapolated data.



## Materials of Manufacture

Filter media	: Polypropylene
Support layer	: Polypropylene
Inner core	: Polypropylene
Outer support	: Polypropylene
End fittings	: Polypropylene
Support ring	: Stainless steel

## Cartridge Dimensions (Nominal)

Diameter	: 70mm (2.8")
Length	: 1 module (short) : 125mm (5")
	1 module : 254mm (10"), 508mm (20")
	2 module : 762mm (30"), 1016mm (40")

## Effective Filtration Area

Up to 0.6m<sup>2</sup> per 250mm module (depending on pore rating).

## Cartridge Treatment

Standard	: Cleaned without further treatment.
Flushed	: Flushed with pyrogen-free water.
Rinsed	: Ultra-clean, pulse flushed to give a system resistivity of 18MΩ.cm

## Gaskets and O-Rings

Ethylene Propylene, FEP encapsulated, Silicone, Viton®, Nitrile or Polypropylene felt.

## Maximum Differential Pressure

### Normal flow direction at :

20°C (68°F)	: 6.0bar (87psi)
80°C (176°F)	: 4.0bar (58psi)
100°C (212°F)	: 3.0bar (44psi)
120°C (248°F)	: 2.0bar (29psi)
125°C (257°F)	: 1.5bar (22psi)

### Reverse flow direction at :

20°C (68°F)	: 2.1bar (30psi)
80°C (176°F)	: 1.0bar (15psi)
100°C (212°F)	: 0.5bar (7psi)

## Operating Temperature

Maximum continuous : 80°C (176°F)

## Sterilisation

*In situ* steam 80 x 30 minute cycles at 135°C (275°F).  
Hot water 200 x 20 minute cycles at 85-90°C (185-194°F).

## Extractables

Minimum total extractables. Please refer to the BorsoPleat-P Validation Guide.

## Integrity Testing

BorsoPleat-P filter cartridges are batch tested for integrity using the Bubble Point Test. Procedural details are available from Van Borselen Filters.




## Clean Water Flow Rates

- Typical clean water flow rate:

A 254mm (10") BorsoPleat-P single cartridge exhibits the flow-ΔP characteristics indicated below, for solutions with a viscosity of 1 centipoise.

- Other solutions:

For solutions with a viscosity of greater than 1 centipoise, multiply the indicated differential pressure by the viscosity in centipoise.



### Range

Suitable for use in **Van Borselen Filters** and as direct replacements for existing cartridges, BorsoPleat-P cartridges can be supplied with end fittings to suit most hardware installations without modification. They are available in single or multiple module units of 5, 10, 20, 30 and 40 inches, and in a choice of removal ratings from 0.5 to 105 micron. BorsoPleat-P Junior versions are also available. Each cartridge is supplied with all necessary seals or O-rings to ensure chemical compatibility.

### Quality Assurance

BorsoPleat-P cartridges are manufactured in an ISO Cleanroom environment by staff fully gowned to minimise any risk of contamination during production. BorsoPleat-P cartridges are batch tested and flushed with pyrogen-free ultra-pure water. As a further safeguard, every cartridge is identified with a batch serial number, allowing users to maintain their own process records.

Registered to ISO 9001, Van Borselen Filters procedures are subject to high standards of quality assurance as demonstrated through its Drug Master File status.

### Material Conformity and Validation

The bio-safety of all materials used in the manufacture of BorsoPleat-P cartridges is assured by FDA approval, USP Class VI and meets or exceeds the latest EC Directives for Food Contact.

A comprehensive validation guide for BorsoPleat-P cartridges is available on request.

### Chemical Compatibility

The BorsoPleat-P materials of construction are compatible with a wide range of chemicals and solvents, however care must be taken to select the appropriate seal material. A comprehensive chemical compatibility guide is available. Since operating conditions vary considerably between applications, verification by the end user is recommended.

### Filter Housings

Please contact a **Van Borselen Filters** representative for further information on our range of filter housings.

#### Van Borselen Filters

Tel: +31 (0)79-3412314

Fax: 079-3412892

PO Box: 3 - 2700 AA Zoetermeer

E-mail: [info@vanborselen.nl](mailto:info@vanborselen.nl)

KVK nr.: Register Haaglanden 27103715

BTW/VAT-Nr.: NL005625452B01

Argonstraat 66, Zoetermeer

[www.vanborselen.nl](http://www.vanborselen.nl)

VAN BORSELEN FILTERS

